



Equipment Disinfection

Product information



Uniclean®



Protect your delicate surface by an active foam

Antistatic foam

Powerful cleaner

Effective in 15 min.

Dries without trace

Broad spectrum

Optimal compatibility



Antistatic



Broad Spectrum



Pleasant Fragrance



Optimal Compatibility



Equipments Disinfection

Uniclean

Foam to clean, disinfect, protect and maintain delicate materials.

Application

- All electronic devices.
- Monitors.
- O.R tables.
- Vinyl surfaces.
- Air conditioners.
- All sensitive surfaces.

Active ingredients

- Propanol.
- Amidon.
- Xylenol.
- Chlorine.

Advantage

- Antistatic.
- Active cleaning foam.
- Powerful cleaner.
- Wipes clean without trace.
- Optimal compatibility with all materials.
- Pleasant fresh lemon scent.
- Fully biodegradable.

Spectrum

Broad-spectrum disinfectant.

- Bactericide.
- Fungicide.
- Virucide(HBV-HCV-HIV).

References

- Bactericidal activity according to the European standard EN 1040.
- Bactericidal activity in 5 minutes, in the presence of interfering substances (dirty conditions) according to the European standard EN 13727.
- Bactericidal activity in 60 minutes, according to the European standard prEN 14561 (carrier test).
- Bactericidal activity in 5 minutes, according to the European standard prEN 14561 (carrier test).
- Tuberculocidal activity against Mycobacterium tuberculosis (BACTEC method). Laboratoire Marcel MERIEUX.
- Fungicidal activity against Candida albicans according to the European standard EN 1275.
- Fungicidal activity against Aspergillus niger according to the European standard EN 14562.
- Efficacy against infectious power of HIV-1 virus (etiologic agent of AIDS) TEXCELL - PARIS
- Activity against PRV, Surrogate of Hepatitis B virus. TEXCELL - PARIS
- Activity against BVDV, Surrogate of Hepatitis C virus. MIKROLAB-Dr J.STEINMANN
- Activity against Herpes virus. MIKROLAB-Dr J.STEINMANN
- Activity against Rotavirus SA 11. MIKROLAB-Dr J.STEINMANN
- Activity against Feline Calicivirus, surrogate of Norovirus. MIKROLAB-Dr J.STEINMANN
- Activity against Influenzavirus H5N1. Retroscreen London.

Method of use



Uniclean is ready for use.

- 1- Pour Uniclean from 1L into the 200ml Foamer bottle.
- 2- Spray the foam onto the surface to be treated.
- 3- Wipe dry with a cloth.

Available in



1Liter With Foamer 200ml.

